

SDK User Guide

Version 2.8.6

InstruStar Electronic Technology

2021-3-21

<http://www.instrustar.com/>

Update

V2.1 (2014.5.13)

- 1, New design

V2.2 (2014.7.10)

- 1, Modified instructions of software trigger description

V2.3 (2014.11.28)

- 1, Modified “ReadVoltageDatas” Time is too long bug

V2.4 (2014.12.5)

- 1, Repair the console class program, unable to detect USB plug problem

V2.5 (2015.7.27)

- 1, Support ISDS2062

V2.6 (2015.8.15)

- 1, Support Roll Mode(Need hardware support)

V2.7 (2016.5.5)

- 1, Support Trigger Sense Div setting (Need hardware support)
- 2, Support Force (Need hardware support)
- 3, Support the setting time of Pulse Width(Need hardware support)
- 4, Support the setting Pre-trigger Percent (Need hardware support)

V2.8 (2017.2.10)

- 1, Support Read voltage datas is out range or not
- 2, Fix 210 Pre-trigger Percent bug
- 3, Add Equipment Id API

V2.8.2 (2019.12.17)

- 1, Support Software control DDS Bias and Zoom

V2.8.6 (2021.3.31)

- 1, Fix restart DLL bug

contents

1. Introduction.....	1
2. Initialization and finish.....	1
3. Equipment Info.....	1
4. Equipment Monitor.....	1
5. Capture Range Set.....	2
6. Sample.....	2
7. Trigger(hardware trigger).....	3
8. AC/DC.....	4
9. Capture Completion Notice.....	4
10. Capture.....	5
11. Data Read.....	5
12. DDS.....	6

1. Introduction

SDK as a virtual oscilloscope equipped with a Windows standard DLL interface, through this interface can directly control a virtual oscilloscope and obtain the data acquired by the oscilloscope. The SDK supports MDSO, MDSO-LA, HDSO, DDSO, ISDS205, ISDS210, ISDS220 and ISDS2062 equipment.

2.Initialization and finish

Call InitDll () to complete the initialization of dynamic library, initialize memory and resources allocated for equipment monitoring and data reading.

int InitDll(void);

Description Dll initialization

Input: -

Output: **Init Status**

Return value 1 Success

0 Failed

Call FinishDll () to complete the end of a dynamic library, freeing memory initialization and related resources in the application.

int FinishDll(void);

Description Dll finished

Input: -

Output: **-Finished Status**

Return value 1 Success

0 Failed

3.Equipment Info

The device ID is a 64-bit integer.

int GetOnlyId0(void);

Description This routines return device id(0-31)

Input: -

Output: - **Device ID(0-31)**

int GetOnlyId1(void);

Description This routines return device id(32-63)

Input: -

Output: - **Device ID(32-63)**

int ResetDevice(void);

Description This routines reset device

Input: -

Output: - **Return value** 1 success

0 failed

4. Equipment Monitor

when the device is detected, the dll have three ways to notify the main program, callback function, set Event and the main program loop detection.

4.1 callback function

When equipment is detect a function "**addcallback**" in the application can be called; when equipment is removed a function "**rmvcallback**" in the application can be called. The DLL has a function pointer which has to be set to these function, using

void SetDevNoticeCallBack(void* ppara, AddCallBack addcallback, RemoveCallBack rmvcallback);

Description This routines sets the callback function of equipment status changed.

Input: **ppara** the parameter of the callback function
 addcallback a pointer to a function with the following prototype:
 void AddCallBack(void * **ppara**)
 rmvcallback a pointer to a function with the following prototype:
 Void RemoveCallBack(void * **ppara**)

Output -

4.2 Event

When equipment is detect, an event "**addevent**" can be set by the DLL; when equipment is removed, an event "**rmvevent**" can be set by the DLL. The user must reset the event when the used them. The DLL has a function pointer which has to be set to these event using

void SetDevNoticeEvent(HANDLE addevent, HANDLE rmvevent);

Description This routines set the event handle, these will be set, when equipment status changed.

Input: **addevent** the event handle
 rmvevent the event handle

Output -

4.3 loop detect

int IsDevAvailable();

Description This routines return the device is available or not.

Input: -

Output **Return value** 1 available
 0 not available

Note: Only need to use one of three ways. Callback and Event functions are asynchronous, more efficient; main program loop detection over a certain time needed to detect whether the device is inserted or removed.

5. Capture Range Set

Device with a programmable gain amplifier, when the signal acquisition time is less than the AD range, the signal amplification gain amplifier to use more AD digits, improving the quality of signal acquisition. Dll will adjusted the range of settings according to the pre-gain amplifier automatically.

int SetOscChannelRange(int channel, int minmv, int maxmv);

Description This routines set the range of input signal.

Input:	channel	the set channel
	0	channel 1
	1	channel 2
	minmv	the minimum voltage of the input signal (mV)
	maxmv	the maximum voltage of the input signal (mV)
Output	Return value	1 Success
		0 Failed

Note: The maximum range of the probe collection X1, the maximum voltage oscilloscope can capture. Like ISDS220 is[-16000mV,16000mV].

Note: In order to achieve better waveform, you need to set the acquisition range, based on the magnitude of the measured waveform. When necessary, you can dynamically change the acquisition range.

6.Sample

int GetOscSupportSampleNum();

Description This routines get the number of samples that the equipment support.

Input: -

Output **Return value** the support sample number

int GetOscSupportSamples(unsigned int* sample, int maxnum);

Description This routines get support samples of equipment.

Input: **sample** the array store the support samples of the equipment
maxnum the length of the array

Output **Return value** the sample number of array stored

int SetOscSample(unsigned int sample);

Description This routines set the sample.

Input: **sample** the set sample

Output **Return value** 0 Failed
other value new sample

7.Trigger(hardware trigger)

This feature requires hardware trigger support. The hardware trigger point is the intermediate data, such as the acquisition of 128K data, trigger point is the 64K point.

Trigger Mode

```
#define TRIGGER_MODE_AUTO 0
#define TRIGGER_MODE_LIANXU 1
```

Trigger Style

```
#define TRIGGER_STYLE_NONE 0x0000 //not trigger
#define TRIGGER_STYLE_RISE_EDGE 0x0001 //Rising edge
#define TRIGGER_STYLE_FALL_EDGE 0x0002 //Falling edge
#define TRIGGER_STYLE_EDGE 0x0004 //Edge
#define TRIGGER_STYLE_P_MORE 0x0008 //Positive Pulse width(>)
```

```

#define TRIGGER_STYLE_P_LESS 0x0010    //Positive Pulse width(>)
#define TRIGGER_STYLE_P      0x0020    //Positive Pulse width(<=)
#define TRIGGER_STYLE_N_MORE 0x0040    //Negative Pulse width(>)
#define TRIGGER_STYLE_N_LESS 0x0080    //Negative Pulse width(<=)
#define TRIGGER_STYLE_N      0x0100    //Negative Pulse width(<=)

```

int IsSupportHardTrigger();

Description This routines get the equipment support hardware trigger or not .

Input: -

Output **Return value** 1 support hardware trigger
0 not support hardware trigger

unsigned int GetTriggerMode();

Description This routines get the trigger mode.

Input: -

Output **Return value** TRIGGER_MODE_AUTO
TRIGGER_MODE_LIANXU

void SetTriggerMode(unsigned int mode);

Description This routines set the trigger mode.

Input: **mode** TRIGGER_MODE_AUTO
TRIGGER_MODE_LIANXU

Output -

unsigned int GetTriggerStyle();

Description This routines get the trigger style.

Input: -

Output **Return value** TRIGGER_STYLE_NONE
TRIGGER_STYLE_RISE_EDGE
TRIGGER_STYLE_FALL_EDGE
TRIGGER_STYLE_EDGE
TRIGGER_STYLE_P_MORE
TRIGGER_STYLE_P_LESS
TRIGGER_STYLE_P
TRIGGER_STYLE_N_MORE
TRIGGER_STYLE_N_LESS
TRIGGER_STYLE_N

void SetTriggerStyle(unsigned int style);

Description This routines set the trigger style.

Input: **style** TRIGGER_STYLE_NONE
TRIGGER_STYLE_RISE_EDGE
TRIGGER_STYLE_FALL_EDGE
TRIGGER_STYLE_EDGE
TRIGGER_STYLE_P_MORE

TRIGGER_STYLE_P_LESS
 TRIGGER_STYLE_P
 TRIGGER_STYLE_N_MORE
 TRIGGER_STYLE_N_LESS
 TRIGGER_STYLE_N

Output -

int GetTriggerPulseWidthNsMin();

Description This routines get the min time of pulse width.

Input: -

Output Return min time value of pulse width(ns)

int GetTriggerPulseWidthNsMax();

Description This routines get the max time of pulse width.

Input: -

Output Return max time value of pulse width(ns)

int GetTriggerPulseWidthDownNs();

Description This routines get the down time of pulse width.

Input: -

Output Return down time value of pulse width(ns)

int GetTriggerPulseWidthUpNs();

Description This routines set the down time of pulse width.

Input: down time value of pulse width(ns)

Output -

void SetTriggerPulseWidthNs(int down_ns, int up_ns);

Description This routines set the up time of pulse width.

Input: up time value of pulse width(ns)

Output -

unsigned int GetTriggerSource();

Description This routines get the trigger source.

Input: -

Output **Return value** 0 :channel 1
1 :channel 2

void SetTriggerSource(unsigned int source);

Description This routines set the trigger source.

Input: **source** 0 :channel 1
1 :channel 2

Output -

int GetTriggerLevel();

Description This routines get the trigger level.

Input: -

Output **Return value** level (mV)

void SetTriggerLevel(int level);

Description This routines set the trigger level.

Input: level (mV)

Output

int IsSupportTriggerSense();

Description This routines get the equipment support trigger sense or not.

Input: -

Return value 1 support
0 not support

int GetTriggerSenseDiv();

Description This routines get the trigger sense.

Input: -

Output **Return value** Sense (0-1 div)

void SetTriggerSenseDiv(int sense);

Description This routines set the trigger sense.

Input: Sense (0-1 div)

Output -

Note: The Range of Trigger Sense is 0.1 Div-1.0 Div. 1 Div =(Sample Range Max Value-Sample Range Min Value)/10.0. For example, Sample Range is [-1000,1000]mV, 1Div=(1000--1000)/10.0 =200mV.

bool IsSupportPreTriggerPercent();

Description This routines get the equipment support Pre-trigger Percent or not .

Input: -

Output **Return value** 1 support
0 not support

int GetPreTriggerPercent();

Description This routines get the Pre-trigger Percent.

Input: -

Output **Return value** Percent (5-95)

void SetPreTriggerPercent(int front);

Description This routines set the Pre-trigger Percent.

Input: Percent (5-95)

Output -

int IsSupportTriggerForce();

Description This routines get the equipment support trigger force or not.

Input: -

Return value 1 support
 0 not support

void TriggerForce();

Description This routines force capture once.

Input: -

Output: -

8.AC/DC

int IsSupportAcDc();

Description This routines get the device support AC/DC switch or not.

Input: -

Output **Return value** 0 :support AC/DC switch
 1 :not support AC/DC switch

void SetAcDc(unsigned int channel, int ac);

Description This routines set the device AC coupling.

Input: channel 0 :channel 1
 1 :channel 2
 ac 1 : set AC coupling
 0 : set DC coupling

Output -

int GetAcDc(unsigned int channel,);

Description This routines get the device AC coupling.

Input: channel 0 :channel 1
 1 :channel 2

Output **Return value** 1 : AC coupling
 0 : DC coupling

9.Capture

Call capture function to begin collecting data, **length** is the length you want to capture, using K Units, such as length = 10, is 10K 10240 points. For sample rate greater than or equal the length of the depth of the collection is stored, take the minimum **length** and depth of storage;For the sampling rate is less than the memory depth, take the minimum **length** and one second data collection length. **force_length** can be forced to cancel the limit of only 1 seconds to be collected.

int Capture(int length, char force_length);

Description This routines set the capture length and start capture.

Input: **length** capture length(KB)

Output	Return value	the real capture length(KB)
--------	---------------------	-----------------------------

Output	memory depth of equipment(KB)
1	100
2	200
3	300
4	400
5	500
6	600
7	700
8	800
9	900
10	1000
11	1100
12	1200
13	1300
14	1400
15	1500
16	1600
17	1700
18	1800
19	1900
20	2000
21	2100
22	2200
23	2300
24	2400
25	2500
26	2600
27	2700
28	2800
29	2900
30	3000
31	3100
32	3200
33	3300
34	3400
35	3500
36	3600
37	3700
38	3800
39	3900
40	4000
41	4100
42	4200
43	4300
44	4400
45	4500
46	4600
47	4700
48	4800
49	4900
50	5000
51	5100
52	5200
53	5300
54	5400
55	5500
56	5600
57	5700
58	5800
59	5900
60	6000
61	6100
62	6200
63	6300
64	6400
65	6500
66	6600
67	6700
68	6800
69	6900
70	7000
71	7100
72	7200
73	7300
74	7400
75	7500
76	7600
77	7700
78	7800
79	7900
80	8000
81	8100
82	8200
83	8300
84	8400
85	8500
86	8600
87	8700
88	8800
89	8900
90	9000
91	9100
92	9200
93	9300
94	9400
95	9500
96	9600
97	9700
98	9800
99	9900
100	10000

Output -

10.3 loop detect

int IsDataReady();

Description This routines return the capture is complete or not.

Input: -

Output **Return value** 1 complete
0 not complete

Note: Only need to use one of three ways. Callback and Event functions are asynchronous, more efficient; main program loop detection over a certain time needed to detect whether the capture is complete or not.

11.Data Read

unsigned int ReadVoltageDatas(char channel, double* buffer,unsigned int length);

Description This routines read the voltage datas. (V)

Input: **channel** **read channel** 0 :channel 1
1 :channel 2
buffer the buffer to store voltage datas
length the buffer length

Output **Return value** the read length

int IsVoltageDatasOutOfRange(char channel);

Description This routines return the voltage datas is out range or not.

Input: **channel** **read channel** 0 :channel 1
1 :channel 2

Output **Return value** 0 :not out range
1 :out range

12.DDS

int IsSupportDDSDevice();

Description This routines get support dds or not

Input: -

Output **Return value** support dds or not

int GetDDSSupportBoxingStyle(int* style);

Description This routines get support wave styles

Input: **style** array to store support wave styles

Output **Return value** if style==NULL return number of support wave styles
else store the styles to array, and return number of wave styles

void SetDDSBoxingStyle(unsigned int boxing);

Description This routines set wave style

Input: **boxing** BX_SINE 0x00 //Sine
BX_SQUARE 0x01 //Square
BX_TRIANGULAR 0x02 //Triangular

BX_UP_SAWTOOTH 0x03 //Up Sawtooth
 BX_DOWN_SAWTOOTH 0x04 //Down Sawtooth

Output: -

void SetDDSPinlv(unsigned int pinlv);

Description This routines set frequency

Input: **pinlv** frequency

Output: -

void SetDDSDutyCycle(int cycle);

Description This routines set duty cycle

Input: **cycle** duty cycle

Output: -

void DDSOutputEnable(int enable);

Description This routines enable dds output or not

Input: **enable** 1 enable
 0 not enable

Output: -

int IsDDSOutputEnable();

Description This routines get dds output enable or not

Input: -

Output **Return value** dds enable or not

int IsDDSSupportSoftwareControlZoomBias();

Description This routines get dds output voltage is support software control

Input: -

Output **Return value** support or not

int GetDDSBiasResistanceRangeMin();

Description This routines get the resistance min value of DDS Bias range.

Input:

Output **Return value** 0 Failed
 other value minimum resistance

int GetDDSBiasResistanceRangeMax();

Description This routines get the resistance max value of DDS Bias range.

Input:

Output **Return value** 0 Failed
 other value maximum resistance

void SetDDSBiasResistance(int Resistance);

Description This routines set the resistance value of DDS Bias.

Input: value resistance

Output

int GetDDSBiasResistance();

Description This routines get the resistance value of DDS Bias.

Input:

Output **Return value 0** Failed
other value resistance

int GetDDSZoomResistanceRangeMin();

Description This routines get the resistance min value of DDS Zoom range.

Input:

Output **Return value 0** Failed
other value minimum resistance

int GetDDSZoomResistanceRangeMax();

Description This routines get the resistance max value of DDS Zoom range.

Input:

Output **Return value 0** Failed
other value maximum resistance

void SetDDSZoomResistance(int Resistance);

Description This routines set the resistance value of DDS Zoom.

Input: value resistance

Output

int GetDDSZoomResistance();

Description This routines get the resistance value of DDS Zoom.

Input:

Output **Return value 0** Failed
other value resistance